





High Performance Liquid Chromatography (HPLC) - Purity Test




Salts and Sugars used as buffering agents are not treated as impurities

Chromatogram Peaks			Number of detected peaks: 3 Analytical methods: UV (214-290 nm), HPLC (Purity), and LC-MS (Identity) Purity: 99.24% Compound: Anastrozole Pubchem CID: 2187 Analysis date: 03/21/26 Calculated Molecular Mass - 293.4 Da
Peak.	RT.	Item.	
1.	2.858	Anastrozole	
2.	6.656		
3.	9.908		

Mass Spectrometric Analysis



Analysis performed by:
 Richard Sylvester MD



LOT-ANAS-00004